


DELTA LONG CANNULA

Venous and arterial long peripheral catheter
with Seldinger technique for multi-day therapies



**LONG CATHETERS, 6-8-10 CM, FOR PERIPHERAL
ARTERIAL AND VENOUS ACCESS WITH DIRECT SELDINGER TECHNIQUE**

AVAILABLE VERSIONS: PV1S AND PV2S



PV1S

- polyurethane radiopaque open tip catheter
- single lumen 3FR or 4FR
- lengths: 6-8-10 cm
- 330PSI ml/sec Power Resistant
- safety echogenic needle 18G or 20G x7 cm
- stainless steel Guidewire 0,018" / 0,032" x30 cm with markers
- Venous and Arterial



PV2S

- polyurethane radiopaque open tip catheter
- single lumen 3FR or 4FR
- lengths: 6-8-10 cm
- 330PSI ml/sec Power Resistant
- safety echogenic needle 18G or 20G x 7 cm
- stainless steel Guidewire 0,018"/0,032" x30 cm with markers
- 325PSI needleless valve
- 10 cm extension line - 300PSI Power Injectable
- sutureless medication and protection patch
- Venous and Arterial

Why choosing DELTA LONG CANNULA?

Double destination of use → Venous and Arterial access

Safety introducer needle → Reduction needle stick injuries and relative possible contamination

PV1S VERSION

Polyurethane catheter (PUR), safety needle, guidewire.

CODE	DESCRIPTION	Q.TY
393063403922	DELTA LONG CANNULA PV1S 3FR X 6CM	20
393083404022	DELTA LONG CANNULA PV1S 3FR X 8CM	20
393103404122	DELTA LONG CANNULA PV1S 3FR X 10CM	20
394063403922	DELTA LONG CANNULA PV1S 4FR X 6CM	20
394083404022	DELTA LONG CANNULA PV1S 4FR X 8CM	20
394103404122	DELTA LONG CANNULA PV1S 4FR X 10CM	20

PV2S VERSION

Polyurethane catheter (PUR) - safety needle, guidewire, extension, needlefree connector, sutureless medication and protection patch.

CODE	DESCRIPTION	Q.TY
393063504222	DELTA LONG CANNULA PV2S 3FR X 6CM	10
393083504322	DELTA LONG CANNULA PV2S 3FR X 8CM	10
393103504422	DELTA LONG CANNULA PV2S 3FR X 10CM	10
394063504222	DELTA LONG CANNULA PV2S 4FR X 6CM	10
394083504322	DELTA LONG CANNULA PV2S 4FR X 8CM	10
394103504422	DELTA LONG CANNULA PV2S 4FR X 10CM	10

Delta Med SpA

via Guido Rossa, 20
46019 Viadana (MN), Italy
Tel.: +39 (0)375 785915 - Fax: +39 (0)375 785201
info@deltamed.it - www.deltamed.it

CERTIFICATES:

CE 0123
EN ISO 9001
EN ISO 13485
MDSAP 13485

FDA Establishment Registration Number: 3006846316